

# Meet EyeBOX<sup>®</sup>, the concussion assessment that is:



## Consistent

Clinically-validated and FDA-cleared.



## Inclusive

Has the ability to assess patients from ages 5 to 67.



## Fast

Delivers a score in under 4 minutes.

## What is the EyeBOX test?

The EyeBOX test is a tool that will help your physician find out if your eye movements are consistent with concussions. Results will be used together with other diagnostics tests to determine the best treatment for you.

## Is EyeBOX right for me?

The EyeBOX is for anyone that has experienced a blow to the head, no matter how minor, in the last week. You do not need a pre-injury baseline for the EyeBOX test to work.

## Are there any risks?

There are no known contraindications or risks to completing the test.



## What should you expect with EyeBOX?

Once you've been properly positioned in the chin rest in front of the EyeBOX device, you will watch a short 4-minute video while the device records subtle eye movements. Once the test is completed, sit back and relax. Your test report will be ready immediately. EyeBOX analyzes the data and produces a BOX score, graphs, and information for your physician to assess.